


As Seen In...

PF PRODUCTS FINISHING

PFOne.com



How to Use XRF Testing Instruments to Accurately Measure Zinc Nickel

Q: What are the XRF testing requirements for plating zinc nickel alloy?

A: X-ray fluorescence (XRF) is a non-destructive method for measuring the quality control of zinc and nickel alloy plating thickness. Accuracy of the XRF measurement can be affected by several factors that can influence the depth of the metal alloy measurement (penetration). Some of these factors are: surface roughness, surface cleanliness, surface preparation, and the XRF instrument used. As a result, it is important that every XRF measurement is taken on a clean, smooth, and flat surface. This can be achieved by using a fine sandpaper or a fine abrasive to smooth the surface. The XRF instrument used should be calibrated to the alloy being measured. This can be done by using a known alloy standard. The XRF instrument should be used in a consistent manner. This means that the same XRF instrument should be used for all measurements. The XRF instrument should be used in a consistent location. This means that the XRF instrument should be used in a consistent location. The XRF instrument should be used in a consistent manner. This means that the same XRF instrument should be used for all measurements. The XRF instrument should be used in a consistent location. This means that the XRF instrument should be used in a consistent location.

PLATING CLINIC

THICKNESS	ALLOY
5 µm	100% Ni/C
10 µm	100% Ni/C
15 µm	100% Ni/C
20 µm	100% Ni/C
25 µm	100% Ni/C
30 µm	100% Ni/C

STANDARD TYPE	THICKNESS	ALLOY
Alloy	5 µm	100% Ni/C
Alloy	10 µm	100% Ni/C
Alloy	20 µm	100% Ni/C
Alloy	30 µm	100% Ni/C
Alloy	40 µm	100% Ni/C
Alloy	50 µm	100% Ni/C
Alloy	60 µm	100% Ni/C

Equipment

So You Don't Have To

RELIANT Finishing Systems

of High-Performance Powder Coating Equipment

reliantfinishing.com info@reliantfinishing.com

40 UNDER 40

Meet the future leaders of the finishing industry.

— page 24

How to Use XRF Testing Instruments to Accurately Measure Zinc Nickel

Q. What are the XRF testing requirements for plating zinc nickel alloy?

A. X-ray fluorescence (XRF) is a common technique used for quality control of zinc and zinc nickel alloy plating deposits. Accurate thickness measurement for both zinc and zinc alloy deposits can be obtained using proper thickness standards.

However, there are several factors that can influence the accuracy of the nickel alloy measurement data obtained by various XRF instruments. As a result, it is common for many XRF units to provide at least partially incorrect data regarding alloy measurement. This can indicate an alloy is too high or too low, depending on the instrument, settings and standards. It can indicate the alloy range is too high when it is not, or that the alloy range is good when it is not.

Measurement Distance: This setting represents the distance that the measurement surface of the part is from the X-ray emitter and fluorescence detector. A higher distance allows for measuring larger parts with recessed low current density (LCD) areas. A unit capable of multiple focal lengths is preferred where possible. Lower distances are preferred for hull cell panels or smaller parts (0-2 cm). Higher distances (3-10 cm) may be needed for measuring recessed LCD areas of larger parts. Lower measurement distances generally allow for the most accurate alloy readings. If the focal length is changed to a different setting, the instrument must be recalibrated. On some multifocal machines, it is possible to maintain separate calibrations for different focal lengths without recalibrating. For this case, the exact focal length should be specified for each calibration.

Measurement Time: Short times allow for quicker analysis but give less consistent results. This is more critical for lower thickness areas. While this can vary depending on the instrument, in general, 30-second measurements are needed for best accuracy in lower thickness areas.

Collimator Size: This represents the size of the measurement beam, indicating the size of the measurement area on the part. A larger collimator size can give a larger, more representative sample, and is generally preferred for QC testing. Typical collimator sizes for zinc nickel are 0.6-1.0 mm. A smaller collimator size can be useful when measuring areas with complicated shapes, or when looking for more detailed info. In this case, bigger is better where possible.

Deposit Thickness: Depending on the above factors, as well as other instrument related details, there will be a minimum thickness required for an instrument that will provide accurate alloy results. This can vary by XRF supplier, and this should be monitored by comparison with test panels provided by your supplier. It may be from 0.1 μm to 5 μm , depending on the instrument and settings.

XRF Calibration Standards: It is typical for zinc nickel

platers to have access to zinc XRF thickness standards. These should be used as part of the calibration set. In general, it is recommended that standards are obtained that cover the alloy range of the specification that must be met.

If you are plating high alloy, you want high alloy standards. If you are plating to meet low alloy specifications, you should have low alloy standards. A minimum of three alloy standards are recommended. Instruments should be recalibrated every four weeks. A daily alloy and thickness calibration check of one standard is also recommended. Follow the XRF supplier recommendations for recalibration frequency. The exact thickness of the zinc standards is not as critical. Use valid zinc calibration standards that are currently in use.

For customers testing low alloy, 6-12 percent nickel specifications, parts typically 5-10 μm , see Table 1.

For testing high alloy, 12-16 percent nickel with large thickness variation, see Table 2.

Hull Cell Testing: Maintaining zinc nickel alloy uniformity is critical to providing consistent corrosion protection. Checking alloy and thickness on hull cell panels, either standard hull cells or long hull cells, is required to monitor the alloy range (low vs. high alloy). Use of the long hull cell allows more detailed monitoring of alloy and thickness over a broad current density range. This is valuable for complicated parts having a large range of current densities.

Having a flat panel allows the tester to measure at the lowest measuring distance possible. It is recommended that a separate calibration be set up that makes use of this, with the measurement distance set to as low as possible, typically 0-1 cm. This will maximize accuracy in low thickness areas. Alloy measurements across the panel can be taken over time.

Some instruments allow for setting up an automated program to take measurements at specific locations, representing certain current densities. Use of this feature is recommended where possible to allow for consistent and efficient measurements.

If all the above factors are not properly considered, XRF testing data can be misleading and can provide inaccurate data regarding nickel alloy, especially in low thickness areas.



MATTHEW STAUFFER
Pavco

Matthew is technical service director at Pavco. Visit pavco.com.

TABLE 1 - Low Alloy: 6-12% Nickel, 5-10 µm Thickness

STANDARD TYPE	THICKNESS	ALLOY
Zinc	5 µm	100% zinc
Zinc	10 µm	100% zinc
Zinc	20 µm	100% zinc
Zinc Nickel	4 µm	5% nickel
Zinc Nickel	8 µm	8% nickel
Zinc Nickel	14 µm	12% nickel

TABLE 2 - High Alloy: 12-16% Nickel, Large Thickness Variation

STANDARD TYPE	THICKNESS	ALLOY
Zinc	5 µm	100% zinc
Zinc	10 µm	100% zinc
Zinc	20 µm	100% zinc
Zinc Nickel	5 µm	10% nickel
Zinc Nickel	10 µm	14% nickel
Zinc Nickel	20 µm	17% nickel

Use of appropriate calibration standards for the applicable alloy specification is critical. After your XRF is optimized per the above recommendations, conduct testing with your supplier to best understand what the minimum thickness

must be in order to provide accurate alloy data. This can vary significantly due to the numerous variables detailed above. Just be sure to ignore any alloy data below this minimum accuracy thickness for your machine. ■■

Our Equipment Works Hard

So You Don't Have To

RELIANT
Finishing Systems

U.S. Manufacturer of High-Performance Powder Coating Equipment

888-770-0021 | reliantfinishingsystems.com | info@reliantfinishingsystems.com

Accurate, non-destructive plating thickness measurement with FISCHER instruments

High precision X-ray fluorescence devices
Exact determination of the thickness and composition of metallic coatings

- Cr/Ni/Cu, ENEPIG and Nickel Phosphate
- Phosphorus content in NiP coatings
- Analysis of metal concentration in electroplating baths

Powerful handheld devices
For precise coating thickness measurement

- Plating coatings such as chromium, copper or zinc in small parts such as screws, nuts and bolts
- Measurement of anodic coatings
- Easy data transfer and analysis software

www.fischer-technology.com | info@fischer-technology.com | 860-683-0781

Fischer

Coating Thickness | Material Analysis | Microhardness | Material Testing